

PRODUCT BROCHURE

OPTICS OBB-A1107

Date: 30.04.2025

ENGLISH



N.A. 0,9 / 0,13 „Swing-out” achromatic condenser (with aperture diaphragm)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Condenser OBB-A1107

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1107
Model-Serie	OBB-AC
Assortment	Optics
Product Type	Microscopes
Product group	Polarising microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761164468
Customs tariff number	90119000

Gross weight 0.80 kg

Net weight 0.30 kg

Shipping method Parcel service

Packaging component - by weight - plastics 500 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Construction

Aperture diaphragm	yes
Condenser Type	Abbe
Condenser aperture	0.9/0.13
Condenser mounting	centerable

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×70 mm
------------------------------	-------------

SPARE PARTS

Model	Description
OBB-A1284	Polarisation unit OBB-A1284