

KERN OBB-A6102UK

Professional illumination units guarantee outstanding, uniform and strong illumination



Article codes, product groups & Status information

Model article-code	OBB-A6102UK
Model-Serie	OZB-IR
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	4045761218475
Planned start date	2019-1HJ
Customs tariff number	85395100
CE Label	Yes

Technical data - Primary

Illumination type	LED
Illumination intensity	4.5W

Technical data - Construction

Color tempertaure	6313 K
Dimmable	Yes
Ring diameter	6.3 cm

Technical data - Power supply

Protection Class	Protection class II - protective insulation
------------------	---

Technical data - Environmental conditions

Storage Temp Min	-5 °C
------------------	-------

Functions

Standard

Stotage Temp Max	40 °C
------------------	-------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	180×110×55 mm
Gross weight	0.50 kg
Net weight	0.40 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	125 g
Packaging component - by weight - plastics	0 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

KERN OBB-A6102UK

Professional illumination units guarantee outstanding, uniform and strong illumination

