

PRODUCT BROCHURE

KERN MBC-A01

Date: 19.09.2025

ENGLISH



Mechanical height rod to fix to the scale, measuring range 10 cm – 80 cm, increments of 1 mm

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Mechanical height rod MBC-A01

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	MBC-A01
Model-Serie	MBC
Assortment	KERN
Product Type	Medical products
Product group	Body size measurement
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761133624
Customs tariff number	84239010
CE Label	yes

Technical data - Primary

Measurement range	80 cm
Readout	1 mm

Technical data - Construction

Dimensions (W×D×H)	900×465×180 mm
--------------------	----------------

Technical data - Measuring system

Measurement range	80 cm
Readout	1 mm
Min Length	10 cm
Device type	mechanical

Units	cm;inch;inch/ft
Tolerance Length	5 mm

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	0 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	180×890×140 mm
Gross weight	1.4 kg
Net weight	0.80 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300.65 g
Packaging component - by weight - plastics	300.65 g
Delivery time	1 d

Technical data - Verification

MPG class	1m
GMDN code	37001

FUNCTIONS

Standard



UNIT



1 DAY

SPARE PARTS

Model	Description
RMBC-9001	Cover RMBC-9001

ENGLISH