

PRODUCT BROCHURE

KERN 304-06

Date: 09.03.2026

ENGLISH



Weight set, E1, 1 g - 1 kg, Knob, stainless steel polished, in wooden case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	304-06
Model-Serie	304
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761002586
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 1 kg
Case wood	313-060-100
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	13
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 g - 1 kg
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×180×110 mm
Gross weight	4.4 kg
Net weight	4.4 kg
Shipping method	Parcel service

Packaging component - 0 g
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 963-219

FUNCTIONS

Standard



Option



SERVICES

Model	Description
963-219	DAkkS Calibration KERN 963-219
962-219R	DAkkS Recalibration KERN 962-219R
921-219	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-219
920-714	Marking (Last digits of serial number) KERN 920-714

SPARE PARTS

Model	Description
313-060-100	Box of Weights KERN 313-060-100
315-242	Tweezers KERN 315-242