

PRODUCT BROCHURE

KERN 354-04

Date: 10.03.2026

ENGLISH



Weight set, M2, 1 g - 200 g, Knob, stainless steel fine turned, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	354-04
Model-Serie	354
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main article
EAN code (Model)	404576 1292307
Predecessor Model	354-44
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 200 g
Case wood	353-440-200
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	11
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 g - 200 g
OIML class	M2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	160×98×58 mm
Gross weight	1.2 kg
Net weight	1.2 kg

Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-617
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-617</u>	DAkkS Calibration KERN 962-617
<u>962-617R</u>	DAkkS Recalibration KERN 962-617R
<u>920-712</u>	Marking (Last digits of serial number) KERN 920-712

SPARE PARTS

Model	Description
<u>354-440-200</u>	Box of Weights KERN 354-440-200
<u>315-242</u>	Tweezers KERN 315-242