

OPTICS ODC 893

Date: 28.01.2025

ENGLISH



USB digital microscope

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

USB digital microscope ODC 893

SPECIFICATION

Article codes, product groups & Status information

Model article-code	ODC 893
Model-Serie	ODC-89
Assortment	Optics
Product Type	Microscopes
Product group	Microscope cameras
Product usage type	Spare part
EAN code (Model)	4045761190542
Planned start date	unbestimmt-1HJ
Customs tariff number	85258900
Discontinued model	yes

Technical data - Primary

Camera type	USB Digital Microscope
Camera resolution	2 MP

Sensor type CMOS

Sensor size 1/3.2"

Technical data - Construction

Color capability	Colour
Interfaces	USB 2.0
Operating system	Win 10; Win 11; Win 7; Win 8
Applicable software	external (OXM 902)
Dimensions	35×35×103 mm

Technical data - Power supply

Power Plug	USB A
Protection Class	Class III - Low voltage
Cable attachment	Permanently connected

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	160×158×85 mm
Gross weight	0.50 kg
Net weight	0.40 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	131.984 g
Packaging component - by weight - plastics	438 mg
Packaging component - by weight - styrofoam	13.432 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Verification

CE	yes
WEEE	yes

ACCESSORIES

Model	Description	Price
ODC-A2404	Calibration Specimen Slide with Micrometer KERN ODC-A2404	40.00 €

RELATED ARTICLES DATA

Model name	Camera type	Camera resolution	Sensor type	Sensor size
ODC-A8506				
ODC-A8507				
ODC 893	USB Digital Microscope	2 MP	CMOS	1/3.2"
ODC-A8201				
ODC-A8301				
ODC-A8501				
ODC-A8502				
ODC-A8503				
ODC-A8504				
ODC-A8505				
ODC-A8701				
ODC-A8801				
ODC-A8901				
ODC-A8902				
ODC-A2401				
ODC-A2402				