

KERN 313-036

Weight set, E2, 1 mg - 100 g, Knob, stainless steel polished, in aluminium case



Article codes, product groups & Status information

Model article-code	313-036
Model-Serie	313
Assortment	Scales
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761494503
Predecessor Model	313-36-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg - 100 g
Case aluminium	313-030-600
Weight set case	In aluminium case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes

Technical data - Measuring system

OIML class	E2
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×125×75 mm
------------------------------	---------------

Functions

Gross weight	0,90 kg
Net weight	0,65 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	250 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-302
-------------------	---------

KERN 313-036

Weight set, E2, 1 mg - 100 g, Knob, stainless steel polished, in aluminium case



Standard



Option

