

# KERN 337-01

Individual weight, OIML Class F2, 1 g, Knob, stainless steel fine turned



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 337-01                         |
| Model-Serie           | 337                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML F2 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761006287                  |
| Planned start date    | 2004-1HJ                       |
| Customs tariff number | 84239090                       |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 1 g         |
| Tolerance (OIML) | 0,3 mg      |
| Case aluminium   | 317-010-600 |
| Case plastic     | 347-030-400 |
| Case wood        | 337-010-200 |

## Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | Stainless steel fine turned |
|----------|-----------------------------|

## Technical data - Measuring system

|                              |        |
|------------------------------|--------|
| Nominal value                | 1 g    |
| OIML class                   | F2     |
| OIML conform                 | Yes    |
| Tolerance (OIML)             | 0,3 mg |
| Density (kg/m <sup>3</sup> ) | 7950   |

## Functions

Standard

Volume (cm<sup>3</sup>) 0,125

## Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 6×6×6 mm       |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Technical data - Verification

|                       |    |
|-----------------------|----|
| Verification possible | No |
|-----------------------|----|

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-431 |
|-------------------|---------|

# KERN 337-01

Individual weight, OIML Class F2, 1 g, Knob, stainless steel fine turned

