

KERN 346-88

Individual weight, OIML Class M1, 20 kg, Block, cast iron lacquered



Article codes, product groups & Status information

Model article-code	346-88
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241657
Predecessor Model	346-88-2019a
Successor Model	346-88-2024e2; 346-88-2025e
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 kg
Tolerance (OIML)	1 g
Case aluminium	346-080-600

Technical data - Construction

Dimensions	230×115×139 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	20 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	1 g

Functions

Density (kg/m ³)	7100
Volume (cm ³)	2816,9014

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	240×130×160 mm
Gross weight	20 kg
Net weight	20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-645
-------------------	---------

KERN 346-88

Individual weight, OIML Class M1, 20 kg, Block, cast iron lacquered



Standard

