

## PRODUCT BROCHURE

---

# KERN 347-066-N00

Date: 28.05.2026

ENGLISH



Individual weight (Newton), 50 cN, Hook weight, stainless steel fine turned

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-066-N00                    |
| Model-Serie           | 347-066-N00                    |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 404576 1474970                 |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Tolerance (OIML) | 3 mg        |
| Case plastic     | 347-070-400 |

### Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | stainless steel fine turned |
|----------|-----------------------------|

### Technical data - Measuring system

|                              |        |
|------------------------------|--------|
| OIML class                   | M1     |
| OIML conform                 | no     |
| Tolerance (OIML)             | 3 mg   |
| Density (kg/m <sup>3</sup> ) | 7950   |
| Volume (cm <sup>3</sup> )    | 5,9524 |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 18×18×36 mm    |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 6 d            |

### Optional Services

|                   |          |
|-------------------|----------|
| DAkkS Calibration | 962-636N |
|-------------------|----------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## ACCESSORIES

| Model                       | Description                         |
|-----------------------------|-------------------------------------|
| <a href="#">347-070-400</a> | Plastic Weight Box KERN 347-070-400 |
| <a href="#">315-242</a>     | Tweezers KERN 315-242               |

## SERVICES

| Model                     | Description                        |
|---------------------------|------------------------------------|
| <a href="#">962-636RN</a> | DAkkS Recalibration KERN 962-636RN |
| <a href="#">962-636N</a>  | DAkkS Calibration KERN 962-636N    |

## RELATED ARTICLES DATA

| Model name                  | Tolerance (OIML) | Case plastic |
|-----------------------------|------------------|--------------|
| <a href="#">347-066-N00</a> | 3 mg             | 347-070-400  |