

## PRODUCT BROCHURE

---

# OPTICS OZB-A5301

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? Holders

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- There are two microscope head holders available for these flexible, modular systems. These brackets are suitable for all stereo microscope stands and universal stands (excluding spring loaded arm), to make focusing possible.
- Holder for stereomicroscope stands with adjustable handwheel rotation
- Diameter of the connector for the microscope head: 76 mm
- Diameter of the connector for the stand: 25 mm

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	OZB-A5301
Model-Serie	OZB-H
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761169609
Customs tariff number	90119000

### Technical data - Construction

Holder diameter	76 mm
-----------------	-------

### Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	290×210×135 mm
Gross weight	1.6 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	239.987 g
Packaging component - by weight - plastics	7.9008 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	81.3124 g
Delivery time	1 d

## FUNCTIONS

Standard

