

LED gooseneck illumination OZB-A4515



Function

1. On/Off-Button
2. Dimmer light intensity



Advantages at a glance

The **OZB-A4515** 6-W-LED gooseneck illuminator allows you to optimally illuminate your sample.



Safety informations

The light from the gooseneck illumination is very bright. Therefore, always switch on the illumination before looking through the eyepieces. Avoid looking directly into the LEDs.



Technical details

Illuminance: 6W-LED
 Colour temperature: 5600-6300K
 Length: 300 mm
 Input voltage: AC 100-240V 50/60HZ



Techn. service +49 7433 9933 777