

PRODUCT BROCHURE

KERN 347-136-N00

Date: 28.05.2026

ENGLISH



Individual weight (Newton), 50 N, Hook weight, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-136-N00
Model-Serie	347-136-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1475038
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	250 mg
Case plastic	347-140-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
OIML conform	no
Tolerance (OIML)	250 mg
Density (kg/m ³)	7950
Volume (cm ³)	595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×130 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services

DAkkS Calibration	962-643N
-------------------	----------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
317-290	Gloves KERN 317-290
347-140-400	Plastic Weight Box KERN 347-140-400

SERVICES

Model	Description
962-643RN	DAkkS Recalibration KERN 962-643RN
962-643N	DAkkS Calibration KERN 962-643N

RELATED ARTICLES DATA

Model name	Tolerance (OIML)	Case plastic
347-136-N00	250 mg	347-140-400