


KERN KDP 300-3

Universal weighing platforms for flexible connection in process chains



Article codes, product groups & Status information

Model article-code	KDP 300-3
Model-Serie	KDP
Assortment	KERN
Product Type	Industrial scales
Product group	Digital weighing systems
Product usage type	Spare part;Main article;Accessories/DL
EAN code (Model)	4045761251885
Predecessor Model	KDP 300-3-2024a
Successor Model	KDP 300-3-????e
Planned start date	2022-1HJ
Customs tariff number	84239010
CE Label	Yes

Technical data - Primary

Weighing capacity	350 g
Readability	1 mg
Indicator cable length	1.5 m
Net weight	1.04 kg
Dimensions weighing platform	165×165×73 mm

Technical data - Construction

Material weighing plate	Plastic, with conductive lacquer
Bubble Level	Yes
Revolving screw feet	Yes

Weighing surface (d)	105 mm
Indicator cable length	1.5 m
Dimensions weighing platform	165×165×73 mm
Dimensions display casing	0×0×0 mm
Dimensions mounted completely	165×166×140 mm
IP protection indicator	No IP protection

Technical data - Power supply

Mains adapter/adapter included	CH;EURO;UK;US
Input Voltage	100 V - 240 V, 50 / 60 Hz

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 5 mg
Readability	1 mg
Recommended adjusting weight	350 g (F1)
Repeatability	2 mg
Warm up time	120 min

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×525×165 mm
------------------------------	----------------

KERN KDP 300-3

Universal weighing platforms for flexible connection in process chains



Gross weight	2 kg
Net weight	1.2 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - plastics	0 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - glass	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	0 g
Packaging component - by weight - metal	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-127
-------------------	---------

Functions

Standard



Option



Factory

