

PRODUCT BROCHURE

KERN 323-076

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 mg - 2 kg, Knob, stainless steel polished, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	323-076
Model-Serie	323
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1223776
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 2 kg
---------------	-------------

Case aluminium 313-070-600

Weight set case in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 2 kg
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	278×180×138 mm
Gross weight	6 kg
Net weight	6 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	320 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-406
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-406</u>	DAkkS Calibration KERN 962-406
<u>962-406R</u>	DAkkS Recalibration KERN 962-406R
<u>921-406</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-406
<u>920-707</u>	Marking (Last digits of serial number) KERN 920-707

SPARE PARTS

Model	Description
<u>313-070-600</u>	Aluminium Protective Box KERN 313-070-600