


PRODUCT BROCHURE

KERN 348-01

Date: 10.03.2026

ENGLISH



Milligram weights, M1, 1 mg, Platelet, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Milligram weights, M1, 1 mg, Platelet, stainless steel polished (OIML)
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	348-01
Model-Serie	348
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761240131
Predecessor Model	348-01-2018a
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg
Tolerance (OIML)	0,2 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	3,5×3,5×3 mm
--------------------	--------------

Material stainless steel polished

Technical data - Measuring system

Nominal value	1 mg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	0,2 mg
Density (kg/m ³)	7960

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-651
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>338-090-200</u>	Wooden Box for Set of Weights KERN 338-090-200
<u>347-009-400</u>	Plastic Weight Box KERN 347-009-400
<u>317-009-600</u>	Aluminium Protective Box KERN 317-009-600
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270

SERVICES

Model	Description
<u>962-651</u>	DAkkS Calibration KERN 962-651
<u>962-651R</u>	DAkkS Recalibration KERN 962-651R