

PRODUCT BROCHURE

OPTICS OZB-A45 15

Date: 30.04.2025

ENGLISH



Double gooseneck illumination for best illumination of your sample

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A45 15
Model-Serie	OZB-IF
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Main article
EAN code (Model)	404576 12 18444
Predecessor Model	OZB-A45 15-2017a
Customs tariff number	94054239

Technical data - Primary

Illumination type	LED
Illumination intensity	6W

Technical data - Construction

Color tempertaure	5600-6300 K
Fibre length	300 mm
Dimensions (W×D×H)	120×370×250 mm

Technical data - Power supply

Protection Class	Protection class I - protective conductor
------------------	--

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	300×250×250 mm
Gross weight	2 kg
Net weight	1.8 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	179.952 g
Packaging component - by weight - plastics	1.25665 g
Packaging component - by weight - styrofoam	70.1211 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	yes
------	-----

FUNCTIONS

Standard



SPARE PARTS

Model	Description
OZB-A4853	Replacement LED OZB-A4853
OZB-A4867	Spare gooseneck with LED KERN OZB-A4867