

KERN 326-36

Individual weight, OIML Class F1, 5 kg, Block, stainless steel polished



Article codes, product groups & Status information

Model article-code	326-36
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761232877
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	25 mg
Case aluminium	346-060-600
Case plastic	Not available
Case wood	Not available

Technical data - Construction

Dimensions	150×68×96 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	5 kg
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	25 mg

Functions

Standard

Density (kg/m³) 7950

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×68×96 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-443
-------------------	---------

KERN 326-36

Individual weight, OIML Class F1, 5 kg, Block, stainless steel polished

