

KERN 317-04

Individual weight, OIML Class E2, 10 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	317-04
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761312524
Predecessor Model	317-4-2.020a; 317-4-2.022i
Planned start date	2022-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	10 g
Tolerance (OIML)	0,06 mg
Case aluminium	317-040-600
Case plastic	317-040-400
Case wood	317-040-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	10 g
OIML class	E2
OIML conform	Yes

Functions

Standard

Tolerance (OIML)	0,06 mg
Density (kg/m ³)	7.950
Volume (cm ³)	1,25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	10×10×18 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-334
-------------------	---------

KERN 317-04

Individual weight, OIML Class E2, 10 g, Knob, stainless steel polished



Option

