

PRODUCT BROCHURE

KERN YKI-01

Date: 16.03.2025

ENGLISH



RS-232 / Ethernet adapter KERN YKI-01

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- Adapter to connect balances, force measuring devices etc. with an RS-232 interface to an IP-based Ethernet network
- The adapter allows direct integration of measuring devices into most company networks, with automatic configuration using DHCP or statical IP allocation
- You can access the current measuring data at any time from within the whole network, without having to use additional cables
- The adapter supports a wide range of configurational options using the enclosed administration software or by the web interface
- The driver delivered with the goods allows you to access the measuring device with existing, non-network-compatible applications using a virtual COM port
- Compatible with KERN BalanceConnection software
- Suitable for wall mounting and DIN track mount or free-standing
- Scope of delivery: RS-232 Ethernet adapter with 1× Port RS-232, mains adapter, CD (driver, software, operating instructions), adhesive feet

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YKI-01
Model-Serie	YKI

Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	4045761146624

Successor Model YKI-10; YKI-11

Customs tariff number 84239010

Technical data - Power supply

Mains adapter/adapter included CH;EURO

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 210×150×72 mm

Gross weight 0.40 kg

Net weight 0.20 kg

Shipping method Parcel service

Packaging component - by weight - cardboard 194 g

Delivery time 1 d

FUNCTIONS

Standard

