

Article codes, product groups & Status information

Model article-code	337-161
Model-Serie	337
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761006652
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	800 mg
Case aluminium	Not available
Case plastic	Not available
Case wood	337-161-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	50 kg
OIML class	F2
OIML conform	Yes
Tolerance (OIML)	800 mg
Density (kg/m ³)	7950

Functions

Standard

Volume (cm³) 6250

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	198×198×249 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	17 d

Optional Services

DAkKS Calibration	962-446
-------------------	---------

KERN 337-161

Individual weight, OIML Class F2, 50 kg, Check weight, stainless steel fine turned

