


# KERN 347-087-N00

Individual weight (Newton), 2 N, Slot weight (Type B), stainless steel fine turned



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-087-N00                    |
| Model-Serie           | 347-7-N00                      |
| Assortment            | Scales                         |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761366725                  |
| Planned start date    | 2024-2HJ                       |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

## Technical data - Primary

|                  |       |
|------------------|-------|
| Tolerance (OIML) | 10 mg |
|------------------|-------|

## Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | Stainless steel fine turned |
|----------|-----------------------------|

## Technical data - Measuring system

|                              |         |
|------------------------------|---------|
| OIML class                   | M1      |
| Density (kg/m <sup>3</sup> ) | 8.000   |
| Volume (cm <sup>3</sup> )    | 23,8095 |

## Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 60×60×14,1 mm  |
| Gross weight                 | 0,20 kg        |
| Net weight                   | 0,20 kg        |
| Shipping method              | Parcel service |

Packaging component - by weight - plastics 0 g

Delivery time 6 d

## Optional Services

DAkkS Calibration 962-638N

## Functions

# KERN 347-087-N00

Individual weight (Newton), 2 N, Slot weight (Type B), stainless steel fine turned



## Standard



## Option

