

## PRODUCT BROCHURE

---

# KERN 324-046

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 g - 200 g, Knob, stainless steel polished, in aluminium case

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 324-046              |
| Model-Serie           | 324                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML F1 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 404576 1223882       |
| Customs tariff number | 84239090             |
| Country of Origin     | DE                   |

### Technical data - Primary

|                 |                   |
|-----------------|-------------------|
| Nominal value   | 1 g - 200 g       |
| Case aluminium  | 314-040-600       |
| Weight set case | in aluminium case |

### Technical data - Construction

|               |                          |
|---------------|--------------------------|
| Material      | stainless steel polished |
| Is weight set | yes                      |
| Type of case  | in aluminium case        |

### Technical data - Measuring system

|               |             |
|---------------|-------------|
| Nominal value | 1 g - 200 g |
| OIML class    | F1          |
| OIML conform  | yes         |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 190×130×80 mm  |
| Gross weight                 | 1.2 kg         |
| Net weight                   | 0.90 kg        |
| Shipping method              | Parcel service |

Packaging component - 200 g  
by weight - plastics

Delivery time 1 d

## Technical data - Verification

Verification possible yes

## Optional Services

DAkkS Calibration 962-417

## FUNCTIONS

Standard



Option



## SERVICES

| Model                           | Description                                          |
|---------------------------------|------------------------------------------------------|
| <a href="#"><u>962-417</u></a>  | DAkkS Calibration KERN 962-417                       |
| <a href="#"><u>962-417R</u></a> | DAkkS Recalibration KERN 962-417R                    |
| <a href="#"><u>921-417</u></a>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-417 |
| <a href="#"><u>920-712</u></a>  | Marking (Last digits of serial number) KERN 920-712  |

## SPARE PARTS

| Model                              | Description                               |
|------------------------------------|-------------------------------------------|
| <a href="#"><u>314-040-600</u></a> | Aluminium Protective Box KERN 314-040-600 |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242                     |