

Personal floor scale with BMI function

Features

- Sturdy steel weighing plate with a non-slip and wear-resistant surface
- Easy and hygienic cleaning
- Secure and non-slip positioning with rubber feet
- Hold function: While weighing patients that are unable to stand still, a mean average weight value is determined
- Simple and convenient 4-key operation
- BMI function to determine underweight/normal weight/surplus weight
- Wall mount for the display device, as standard, only for models without a stand
- Battery- or mains-powered
- Protective working cover included with delivery


Technical data


- LCD display, digit height 25 mm
- Dimensions of display device W×D×H 210×110×50 mm
- KERN MPB: Cable length of display device approx. 1,2 m
- Overall dimensions W×D×H
KERN MPB: 315×300×60 mm
KERN MPB-P: 315×450×1000 mm (incl. stand)
- Battery operation possible, 6×1.5 V AA not included, operating time up to 170 h
- Mains adapter external, standard


Accessories


- **1** Mechanical height rod, measuring range: 60 cm – 205 cm, can be screwed onto the column from the front or the back (MPB-P) or wall-mounted, approved as a medical device according to 93/42/EEC or regulation (EU) 2017/745, KERN MSF 200
- Stand to elevate display device, height of stand approx. 950 mm, KERN MWS-A01
- Alcohol-free cleaning cloths for disinfectant wiping, quick acting, based on modern quaternary ammonium compounds and effective against papova viruses. Particularly gentle on materials. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 100 pcs., KERN MYC-01

STANDARD

 MOVE

 BATT

 MULTI

 1 DAY

Model	Weighing range	Readout	Weighing plate	Net weight
	[Max] kg	[d] kg	W×D×H mm	approx. kg
KERN				
MPB 300K100	300	0,1	315×300×60	4,6
MPB 300K100P	300	0,1	315×300×55	7