

PRODUCT BROCHURE

OPTICS OBB-A1240

Date: 16.03.2025

ENGLISH



Infinity planachromatic objective lens 100 x / 1,25 W.D., (0,19 mm) (oil) (sprung)

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Infinity planachromatic objective lens 100 x / 1,25 W.D., (0,19 mm) (oil) (sprung) OBB-A1240

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OBB-A1240
Model-Serie	OBB-AO
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part; Accessories/Services
EAN code (Model)	4045761167728
Customs tariff number	90119000

FUNCTIONS

Standard



ENGLISH

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	65×40×40 mm
Gross weight	0.60 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	300 g
Delivery time	1 d