

KERN 347-076

Individual weight, OIML Class M1, 100 g, Hook weight, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-076
Model-Serie	347-076
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256217
Predecessor Model	347-476
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	5 mg
Case plastic	347-090-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	100 g
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	5 mg
Density (kg/m ³)	7.950

Functions

Standard

Volume (cm³) 11,9048

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	22×22×54 mm
Gross weight	0,10 kg
Net weight	0,10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-637
-------------------	---------

KERN 347-076

Individual weight, OIML Class M1, 100 g, Hook weight, stainless steel fine turned



Option

