

PRODUCT BROCHURE

KERN 356-77

Date: 23.03.2025

ENGLISH



Individual weight, OIML Class M2, 10 kg, Block Eco, cast iron lacquered

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	356-77
Model-Serie	356
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241794
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	1,6 g
Case aluminium	346-070-600

Technical data - Construction

Dimensions	190×95×109 mm
Material	cast iron lacquered

Technical data - Measuring system

Nominal value	10 kg
OIML class	M2
OIML conform	yes
Tolerance (OIML)	1,6 g
Density (kg/m ³)	7100
Volume (cm ³)	1408,4507

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	190×95×109 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-644
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
346-070-600	Aluminium Protective Box KERN 346-070-600

SERVICES

Model	Description
962-644	DAkkS Calibration KERN 962-644
962-644R	DAkkS Recalibration KERN 962-644R
952-644	Verification KERN 952-644
920-30	Marking with stroke number KERN 920-30
920-301	Marking with stroke number (Last digits of serial number) KERN 920-301

ENGLISH

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
356-86	5 kg	800 mg	346-060-600
356-76	5 kg	800 mg	346-060-600
356-77	10 kg	1,6 g	346-070-600
356-87	10 kg	1,6 g	346-070-600
356-88	20 kg	3 g	346-080-600
356-78	20 kg	3 g	346-080-600
356-79	50 kg	8 g	-
356-89	50 kg	8 g	-