

KERN 347-147

Individual weight, OIML Class M1, 10 kg, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-147
Model-Serie	347-147
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366879
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	500 mg

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	500 mg
Density (kg/m ³)	8.000
Volume (cm ³)	1.190,4762

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	178×178×62 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-644
-------------------	---------

Functions

KERN 347-147

Individual weight, OIML Class M1, 10 kg, Slot weight (Type B), stainless steel fine turned



Standard



Option

