

PRODUCT BROCHURE

KERN 344-086

Date: 10.03.2026

ENGLISH



Weight set, M1, 1 g - 5 kg, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	344-086
Model-Serie	344
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224827
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 5 kg
Case aluminium	314-080-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	16
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 5 kg
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	460×260×200 mm
Gross weight	18 kg
Net weight	13 kg

Shipping method	Parcel service
Packaging component - by weight - plastics	5.4 kg
Delivery time	1 d

Optional Services

DAkkS Calibration	962-621
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-621	DAkkS Calibration KERN 962-621
962-621R	DAkkS Recalibration KERN 962-621R
920-716	Marking (Last digits of serial number) KERN 920-716

SPARE PARTS

Model	Description
313-080-600	Aluminium Protective Box KERN 313-080-600
315-242	Tweezers KERN 315-242