

EU Declaration of Conformity | EU-Konformitätserklärung

EN Wir erklären hiermit unter alleiniger Verantwortung, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Richtlinien übereinstimmt. Für die Waage liegt eine EG-Bauartzulassung vor. Das Jahr der EG-Eichung ist neben der CE Konformitätskennzeichnung aufgeführt. Solche Waagen sind ab Werk EG-geeicht und tragen die Kennzeichnung „M“ auf dem Gerät selbst und auf der Verpackung. Sie dürfen im eichpflichtigen Verkehr verwendet werden. Die Fallbeschleunigung gibt bei EG-geeichten Waagen an, für welchen Aufstellungsort die Waage EG-geeicht ist.

DE We hereby declare and assume sole responsibility for the declaration that the product complies with the standards hereinafter. The scale was issued with an EC Type Approval Certificate. The year of EC verification can be found next to the CE mark of conformity. Such scales are EC verified ex factory and bear the mark “M”, on the device itself as well as on the packaging. They can be used for weighing with obligation to verify. For EC verified scales, gravity acceleration indicates for which place of installation the scale is EC verified.

Model | Modell |

Type | Typ |


DLB 160-3A
DLB 160-3A110V

TDLG 160-3-A
TDLG 160-3-A110V

CE mark applied CE Kennzeichnung	EU directive EU-Richtlinie	Standards Normen
	2011/65/EU (RoHS)	EN 63000:2018
	2014/30/EU (EMC)	EN 61326-1:2013
	2014/35/EU (LVD)	EN 61010-1:2010 EN 61010-2-010:2019

Date | Datum | : 26.10.2021

Place of issue: 72336 Balingen,
Ort der Ausstellung: Germany



Albert Sauter
KERN & Sohn GmbH

Signature: Managing director
Signatur: Geschäftsführer