

KERN EMB 6000-1

Entry level laboratory balance with tremendous weighing performance



Article codes, product groups & Status information

Model article-code	EMB 6000-1
Model-Serie	EMB
Assortment	KERN
Product Type	Laboratory balances
Product group	School balances
Product usage type	Main article
EAN code (Model)	4045761111387
Planned start date	2010-1HJ
Customs tariff number	84238129
CE Label	Yes

Technical data - Primary

Weighing capacity	6000 g
Readability	0.1 g
Repeatability	0.1 g
Linearity	± 0.3 g
Weighing surface (d)	150 mm
Recommended adjusting weight	5 kg (M1)

Technical data - Construction

Dimensions of housing	170×244×39 mm
Material weighing plate	Plastic
Casing material	Plastic
Weighing surface (d)	150 mm

Technical data - Display

Display height	15 mm
----------------	-------

Backlight	No
-----------	----

Technical data - Power supply

Battery	9 V Block
Operating time	25 h

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 0.3 g
Readability	0.1 g
Recommended adjusting weight	5 kg (M1)
Repeatability	0.1 g
Stabilisation time under laboratory conditions	3 s
Tare range	6000 g
Warm up time	120 min
Weighing system	Strain gauge
Weighing units	Dwt;g;oz;ozt

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	200×265×66 mm
------------------------------	---------------

KERN EMB 6000-1

Entry level laboratory balance with tremendous weighing performance



Gross weight 0.70 kg

Net weight 0.55 kg

Shipping method Parcel service

Packaging component - by weight - paper 5 g

Packaging component - by weight - cardboard 720 g

Packaging component - by weight - plastics 5 g

Packaging component - by weight - styrofoam 98 g

Packaging component - by weight - foam 23 g

Delivery time 1 d

Optional Services

DAkkS Calibration 963-128

Functions

Standard



Option

