

## PRODUCT BROCHURE

---

# OPTICS OBB-A1580

Date: 28.03.2025

ENGLISH



Trinocular lens barrel

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## DESCRIPTION

Trinocular lens barrel OBB-A1580

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	OBB-A1580
Model-Serie	OBB-AT
Assortment	Optics
Product Type	Microscopes
Product group	Transmitted light microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761246485
Predecessor Model	OBB-A1346
Customs tariff number	90119000

### Technical data - Primary

Tube	Trinocular
------	------------

### Technical data - Construction

Tube type	Butterfly
Tube angle	30°
Tube rotation	yes
Eye distance (max)	75 mm
Eye distance (min)	50 mm
Adjustable dioptre	one-sided

### Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	265×185×165 mm
Gross weight	2 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	600 g
Delivery time	1 d

## FUNCTIONS

Standard
----------

