

KERN 367-12

Individual weight, OIML Class M3, 2 kg, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	367-12
Model-Serie	367-12
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761257306
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	2 kg
Tolerance (OIML)	1 g
Case aluminium	317-120-600
Case plastic	347-120-400
Case wood	337-120-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	2 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	1 g

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	238,0952

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×100 mm
Gross weight	2,0 kg
Net weight	2,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkkS Calibration	962-642
-------------------	---------

KERN 367-12

Individual weight, OIML Class M3, 2 kg, Knob, stainless steel fine turned



Option

