

PRODUCT BROCHURE

KERN 334-086

Date: 10.03.2026

ENGLISH



Weight set, F2, 1 g - 5 kg, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	334-086
Model-Serie	334
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	4045761224421
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 5 kg
Case aluminium	314-080-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 5 kg
OIML class	F2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	460×260×200 mm
Gross weight	18 kg
Net weight	13 kg
Shipping method	Parcel service

Packaging component - 5.4 kg
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 962-421

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-421	DAkkS Calibration KERN 962-421
962-421R	DAkkS Recalibration KERN 962-421R
921-421	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-421
920-716	Marking (Last digits of serial number) KERN 920-716

SPARE PARTS

Model	Description
313-080-600	Aluminium Protective Box KERN 313-080-600
315-242	Tweezers KERN 315-242