

PRODUCT BROCHURE

KERN 354-096

Date: 10.03.2026

ENGLISH



Weight set, M2, 1 g - 10 kg, Knob, stainless steel fine turned, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	354-096
Model-Serie	354
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main article
EAN code (Model)	404576 1292468
Predecessor Model	354-496
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 10 kg
Case aluminium	314-090-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 10 kg
OIML class	M2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	610×260×235 mm
Gross weight	28 kg
Net weight	24 kg
Shipping method	Parcel service

Packaging component - 4.1 kg
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 962-622

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-622	DAkkS Calibration KERN 962-622
962-622R	DAkkS Recalibration KERN 962-622R
920-717	Marking (Last digits of serial number) KERN 920-717

SPARE PARTS

Model	Description
313-090-600	Box of Weights KERN 313-090-600
315-242	Tweezers KERN 315-242