


PRODUCT BROCHURE

KERN 346-055-004

Date: 17.04.2025

ENGLISH



Weight container for rectangular weights, calibrated to OIML class M1

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Preconfigured weight container for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container is adjusted to OIML accuracy class M1. Other OIML accuracy classes are also available, please ask
- Weight of the weight container, OIML class M1: 50 kg
- Possible equipment, rectangular weights, OIML class M1: 4×50 kg
- Maximum total weight (weight container incl. weights): 250 kg
- Tolerance $\pm 2,5$ g

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	346-055-004
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1225084
Customs tariff number	84239090

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	2,5 g

Technical data - Construction

Dimensions	1000×500×500 mm
------------	-----------------

Material stainless steel

Technical data - Measuring system

Nominal value	50 kg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	2,5 mg

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	810×380×390 mm
Gross weight	55 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	5 kg
Delivery time	2 d

Optional Services

DAkkS Calibration 962-646

FUNCTIONS

Standard



SERVICES

Model	Description
962-646	DAkkS Calibration KERN 962-646
962-646R	DAkkS Recalibration KERN 962-646R

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-06	5 kg	250 mg	346-060-600
346-07	10 kg	500 mg	346-070-600
346-08	20 kg	1 g	346-080-600
346-09	50 kg	2,5 g	346-090-600
346-022-005	20 kg	1 g	-
346-042-008	40 kg	1,5 g	-
346-052-010	50 kg	2,5 g	-
346-055-004	50 kg	2,5 g	-
346-055-009	50 kg	2,5 g	-
346-065-009	60 kg	3 g	-
346-86	5 kg	250 mg	346-060-600
346-87	10 kg	500 mg	346-070-600
346-77	10 kg	500 mg	346-070-600
346-76	5 kg	250 mg	346-060-600
346-79-2023a	50 kg	2,5 g	-
346-88	20 kg	1 g	346-080-600
346-78	20 kg	1 g	346-080-600
346-89	50 kg	2,5 g	-

Model name	Nominal value	Tolerance (OIML)	Case aluminium
346-81	100 kg	5 g	-
346-82	200 kg	10 g	-
346-83	500 kg	25 g	-
346-84	1000 kg	50 g	-
346-79	50 kg	2,5 g	346-090-600