

PRODUCT BROCHURE

KERN 325-062

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 mg - 1 kg, Eco Shape, stainless steel polished, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	325-062
Model-Serie	325
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 122426 1
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 1 kg
---------------	-------------

Case wood	313-060-100
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 mg - 1 kg
OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×180×110 mm
Gross weight	2.4 kg
Net weight	2.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-405
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-405	DAkkS Calibration KERN 962-405
962-405R	DAkkS Recalibration KERN 962-405R
921-405	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-405
920-706	Marking (Last digits of serial number) KERN 920-706

SPARE PARTS

Model	Description
313-060-100	Box of Weights KERN 313-060-100
315-242	Tweezers KERN 315-242

Model	Description
<u>318-270</u>	Dusting Brush KERN 318-270