

PRODUCT BROCHURE

KERN 315-274

Date: 10.03.2026

ENGLISH



Weight grip, plastic coated

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Weight grip, plastic coated
- Suitable for knob shaped weights with weight value 5 kg of classes E1 - M3
- Not appropriate for cast iron weights

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	315-274
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	Cases, boxes & accessories
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1222960
Customs tariff number	84239090
Country of Origin	CN

Technical data - Packing & Shipping

Packaging component - by weight - plastics	100 g
Delivery time	1 d

FUNCTIONS

Standard



SERVICES

Model	Description
970-010	Warranty extension KERN (+3 years) 970-010