

PRODUCT BROCHURE

KERN
SET-51

Date: 18.03.2025

ENGLISH

Software EasyTouch SET-51 **Audit Log**

EasyTouch Audit Logs – for critical weighing processes

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Note: Prerequisite for this set is the basic program SET-01 Base
- The SET-51 Audit Logs function logs all actions and activities in EasyTouch to guarantee full monitoring
- Audit Trail: All actions will be stored in the system in a tamper-proof memory and are audit trail compliant, to fulfil the highest security standards. This is vital in critical sectors such as the pharmaceutical industry, the manufacture of food-related products and the automotive and aviation industries, where comprehensive documentation of weighing processes is essential
- Login and Logout: Each login and logout in EasyTouch will be logged with the user name and time stamp
- User administration: All actions relating to the creation, modification or deletion of users will be reliably logged in the Audit Log
- User roles: User role administration will be documented, including creation, modification or deletion of roles
- Master data management: Changes to master data such as creation, modification or deletion will also be recorded
- Logging of measuring data: All measurement and weighing results from all active EasyTouch functions (here from Batch and Statistics) are also recorded in audit logs
- Software settings: Modifications to the general settings, including the Admin area, will be noted in the Audit Log

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	SET-51
Model-Serie	ETS
Assortment	KERN
Product Type	Software
Product group	EasyTouch
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1486355
Customs tariff number	85234910

Technical data - Packing & Shipping

Delivery time	1 d
---------------	-----

FUNCTIONS

Standard



Option

