

KERN OZP 152

Individuality, variety and flexible working through our modular construction system ?
Stereo microscope heads



Article codes, product groups & Status information

Model article-code	OZP 152
Model-Serie	OZB-M
Assortment	Microscopes & Optical Metrology
Segment	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Main article
EAN code (Model)	4045761504509
Predecessor Model	OZM 547; OZP 552
Planned start date	2027K -1HJ
Customs tariff number	90119000
Country of Origin	CN

Technical data - Primary

Optical System	Greenough
Tube	Trinocular
Objectives	0,65x - 5,5x

Technical data - Construction

System	Upright
Tube angle	45°
Tube rotation	Yes
Eye distance (max)	76 mm
Eye distance (min)	54 mm
Adjustable dioptrre	Both-sided
Visual field min	4,2 mm

Functions

Standard

Visual field max	35,4 mm
Working distance	105 mm
Sample height max	90 mm

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	300×250×255 mm
Gross weight	2,6 kg
Net weight	1,8 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	383,2 g
Packaging component - by weight - plastics	12,8 g
Packaging component - by weight - styrofoam	404 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

KERN OZP 152

Individuality, variety and flexible working through our modular construction system ?
Stereo microscope heads

