

**PRODUCT BROCHURE**

---

**KERN  
3 17-06-2020a**

Date: 15.03.2025

**ENGLISH**

Individual weight, OIML Class E2, 50 g, Knob, stainless steel polished

---

**CONTACT****KERN & SOHN GmbH**

Ziegelei 1

72336 Balingen

Germany

Phone : +49 7433 9933-0

Fax : +49 7433 9933-149

Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)Web : [www.kern-sohn.com](http://www.kern-sohn.com)

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	317-06-2020a
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Trade parts list element
EAN code (Model)	4045761003620
Successor Model	317-06; 317-06-2022i
Customs tariff number	84239090

### Technical data - Primary

Nominal value	50 g
Tolerance (OIML)	0,1 mg
Case aluminium	317-060-600
Case plastic	317-060-400
Case wood	317-060-100

### Technical data - Construction

Material	stainless steel polished
----------	--------------------------

### Technical data - Measuring system

OIML class	E2
OIML conform	yes
Density (kg/m <sup>3</sup> )	7950
Volume (cm <sup>3</sup> )	6,25

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	18×18×28 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

### Technical data - Verification

Verification possible	yes
-----------------------	-----

### Optional Services

DAkkS Calibration	962-336
-------------------	---------

### FUNCTIONS

Standard
----------



## ACCESSORIES

Model	Description
317-060-100	Wooden Box for Set of Weights KERN 317-060-100
317-060-400	Plastic Weight Box KERN 317-060-400
317-060-600	Aluminium Protective Box KERN 317-060-600
315-242	Tweezers KERN 315-242

## SERVICES

Model	Description
952-336	Verification KERN 952-336
962-336	DAkkS Calibration KERN 962-336
962-336R	DAkkS Recalibration KERN 962-336R
920-33	Marking by etching KERN 920-33
921-336	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-336
920-331	Marking by etching (last digits of the serial number) KERN 920-331