

PRODUCT BROCHURE

KERN 312-046

Date: 10.03.2026

ENGLISH



Weight set, E2, 1 g - 200 g, Compact form, stainless steel polished, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	312-046
Model-Serie	312
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761223509
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 200 g
Case aluminium	314-040-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes

Technical data - Measuring system

OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	190×130×80 mm
Gross weight	1.2 kg
Net weight	0.90 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	200 g
Delivery time	1 d

Technical data - Verification

Verification possible no

Optional Services

DAkkS Calibration 962-317

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-317</u>	DAkkS Calibration KERN 962-317
<u>962-317R</u>	DAkkS Recalibration KERN 962-317R
<u>921-317</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-317
<u>920-712</u>	Marking (Last digits of serial number) KERN 920-712

SPARE PARTS

Model	Description
<u>314-040-600</u>	Aluminium Protective Box KERN 314-040-600
<u>315-242</u>	Tweezers KERN 315-242