

KERN 327-05

Individual weight, OIML Class F1, 20 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	327-05
Model-Serie	327
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761004726
Successor Model	327-05-2025e
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	20 g
Tolerance (OIML)	0,25 mg
Case aluminium	317-050-600
Case plastic	347-050-400
Case wood	317-050-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	20 g
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	0,25 mg

Functions

Standard

Density (kg/m³) 7950

Volume (cm³) 2,5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	13×13×23 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible Yes

Optional Services

DAkkS Calibration 962-435

KERN 327-05

Individual weight, OIML Class F1, 20 g, Knob, stainless steel polished

