

## PRODUCT BROCHURE

---

# OPTICS OZB-A6303

Date: 21.10.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? ECO-Universal stands

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

Stereomicroscope Stand (Universal) OZB-A6303

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                 |
|-----------------------|---------------------------------|
| Model article-code    | OZB-A6303                       |
| Model-Serie           | OZB-UE                          |
| Assortment            | Optics                          |
| Product Type          | Microscopes                     |
| Product group         | Stereo microscopes              |
| Product usage type    | Spare part;Accessories/Services |
| EAN code (Model)      | 404576 123502 1                 |
| Customs tariff number | 90 11 9000                      |

### Technical data - Construction

|            |                                                                      |
|------------|----------------------------------------------------------------------|
| Stand type | universal (spring loaded arm with clamp/holder/coarse focusing knob) |
|------------|----------------------------------------------------------------------|

Dimensions (W×D×H) 770×85×300 mm

### Technical data - Environmental conditions

|                  |       |
|------------------|-------|
| Storage Temp Min | -5 °C |
| Storage Temp Max | 40 °C |

### Technical data - Packing & Shipping

|                                             |                |
|---------------------------------------------|----------------|
| Dimensions packaging (W×D×H)                | 580×420×370 mm |
| Gross weight                                | 9 kg           |
| Net weight                                  | 7 kg           |
| Shipping method                             | Parcel service |
| Packaging component - by weight - cardboard | 750 g          |
| Packaging component - by weight - styrofoam | 750 g          |
| Delivery time                               | 1 d            |

## FUNCTIONS

Standard



## SPARE PARTS

| Model     | Description                          |
|-----------|--------------------------------------|
| OZB-A6390 | Holder for spring link arm OZB-A6390 |
| OZB-A6392 | Spare clamping lever OZB-A6392       |