

OPTICS OZB-A458 1

Date: 29.01.2025

ENGLISH



Power supply for range for inverse microscope

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

Euro power supply unit OZB-A4581

SPECIFICATION

Article codes, product groups & Status information

Model article-code	OZB-A4581
Model-Serie	OZB-AZ
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Spare part
EAN code (Model)	4045761217751
Planned start date	unbestimmt-1HJ
Customs tariff number	85044095

Technical data - Power supply

Power supply kind	Power adapter
Input Voltage	100 V - 240 V, 50 / 60 Hz

Protection Class	Protection class II - protective insulation
------------------	---

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	200×100×50 mm
Gross weight	2 kg
Net weight	1.00 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	1 kg
Delivery time	1 d

Technical data - Verification

CE	yes
WEEE	yes

RELATED ARTICLES DATA

Model name	Illumination type	Illumination intensity	Weighing capacity
OZB-A4554			
OZB-A4590			
OZB-A4595			
OZB-A4596			
OZB-A4536			
OZB-A4581			
OZB-A4582			
OZB-A7906			
OZB-A7907			
OZB-A7908			
OZB-A4534	LED		
OZB-A4535			
OZB-A4591			
OZB-A7901			
OZB-A7902			
OZB-A7903			
OZB-A7904			
OZB-A7905			

ENGLISH