



Display device: Plastics, protection against dust and water splashes, integrated power supply

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

Display device KMF-TM, LCD display, display digit height 25 mm, weighing units kg

SPECIFICATION

Article codes, product groups & Status information

Model article-code	KMF-TM
Model-Serie	KMF
Assortment	KERN
Product Type	Medical products
Product group	Personal scales
Product usage type	Spare part
EAN code (Model)	4045761150492
Planned start date	2014-1HJ
Customs tariff number	84239010
Discontinued model	yes

Technical data - Display

Display height	25 mm
----------------	-------

Technical data - Power supply

Charging time	12 h
Input Voltage	220 V - 240 V AC, 50 Hz
Operating time (Backlight off)	22 h
Rechargeable battery	Rchrg battery optional - intern

Technical data - Measuring system

Calibration method	Adjusting with external weight
Weighing units	kg

Technical data - Environmental conditions

Maximum humidity	80 %
Maximum operating temperature	40 °C

Minimum ambient temperature	-10 °C
-----------------------------	--------

Minimum humidity	0 %
------------------	-----

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	370×200×95 mm
------------------------------	---------------

Gross weight	0.90 kg
--------------	---------

Net weight	0.70 kg
------------	---------

Shipping method	Parcel service
-----------------	----------------

Packaging component - by weight - cardboard	100 g
---	-------

Packaging component - by weight - plastics	100 g
--	-------

Delivery time	1 d
---------------	-----

Technical data - Verification

Verification approval conforming to 2014/31/EU	yes
--	-----

Verification class	III
--------------------	-----

ACCESSORIES

Model	Description	Price
911-013	Matrix needle printer KERN 911-013	480.00 €

SPARE PARTS

Model	Description	Price
RMCC-5001	Power socket RMCC-5001	5.00 €
RMCC-9002	Display protective plate RMCC-9002	7.00 €
MBC-A09	Holder for batteries KERN MBC-A09	15.00 €
RMPE-8004	Battery compartment cover RMPE-8004	5.00 €

RELATED ARTICLES DATA

Model name	Weighing capacity	Readability	Verification value (e)	Minimum piece weight
RMPEN-4012				
RTMPN-9002				
RMPC-4003				
RMPEN-2002				
RTMPN-9001				
RMPC-9002				
RMPE-8006				
RMPE-8015				
RMPE-2001				
RMPD-2003				
RMPEN-2001				
RMPEN-4011				
RMPCN-4011				
RMPE-9017				
RMPE-8005				
RMPS-8010				
RMPI-8001				
RMPD-2002				
RMPD-8002				
RMPE-9018				
RMPE-6011				

RMPDN-40 10

RMPEN-40 10

RMPCN-40 10

RMPEN-800 1

RMPD-200 1

RMPC-400 1

RMPC-100 2

RMPC-800 2

RMPS-40 10

RMPT-40 10

RMPD-E-100 1 300 kg

RMPD-E-400 1

RMPD-E-200 1

RMPC-100 1

RMPS-400 3

RMPS-800 4

RMPE-90 1 2

RMPD-E-800 1

RMPS-800 7

RMPT-400 1

RMPT-800 5

RMPS-800 8

RMPS-800 9

RMPB-400 1

RMPS-200 3

RMPS-300 1

RMPS-400 1

RMPS-800 3

RMPS-800 6

RMXS-500 1

RMPS-600 1

RMPS-100 1

RMPE-90 1 1

KMF-TM

RMPE-900 8

RMPE-9009

RMPE-9005

RMPD-4001

RMPE-9006

RMPE-8004

RMPE-9007

RMPE-6001

RMPE-9003

RMPE-9004

RMPD-9001

RMPD-9002

RMPC-4002

RMPE-8003
