

KERN ORA 80BE

Refractive index measurement for laboratories and the industry



Article codes, product groups & Status information

Model article-code	ORA 80BE
Model-Serie	ORA-E
Assortment	Microscopes & Optical Metrology
Product Type	Refractometers
Product group	Analogue refractometers
Product usage type	Main article
EAN code (Model)	4045761160149
Planned start date	2015-1HJ
Customs tariff number	90275000
Country of Origin	CN

Technical data - Primary

Scales	Brix
Scale Range	0 % - 50 % 50 % - 80 %

Technical data - Construction

Type	Handheld refractometer
Display scale	Analogue
Dimensions (W×D×H)	180×40×40 mm

Technical data - Measuring system

Measuring method	Light refraction
Measuring Temp min	10 °C
Measuring Temp max	30 °C
Readout	0,5 % 0,5 %

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	215×95×65 mm
Gross weight	0,85 kg
Net weight	0,80 kg
Shipping method	Parcel service
Packaging component - by weight - paper	842,4 mg
Packaging component - by weight - cardboard	50,0175 g
Packaging component - by weight - plastics	1,84275 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Optional Services

ISO Calibration	961-290
-----------------	---------

Functions

KERN ORA 80BE

Refractive index measurement for laboratories and the industry



Standard



Option

