

KERN 347-057-N00

Individual weight (Newton), 20 cN, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-057-N00
Model-Serie	347-7-N00
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366664
Planned start date	2024-2HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	2,5 mg
------------------	--------

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
Density (kg/m ³)	8.000
Volume (cm ³)	2,381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	38×38×4 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 6 d

Optional Services

DAkkS Calibration 962-635N

Functions

KERN 347-057-N00

Individual weight (Newton), 20 cN, Slot weight (Type B), stainless steel fine turned



Standard



Option

