

NEW



OBN 142C861



OBN 142T251

Fluorescence digital set: high-resolution image quality for complex analyses

Features

- The KERN OBN 142 fluorescence microscope is now also available as a complete digital solution in a practical set with camera and matching C-mount adapter. It combines high-quality optics with a stable design and 3 W LED transmitted illumination and 5 W LED epifluorescence incident illumination for excellent image quality
- The sets include either a KERN ODC 861 20 MP fluorescence camera or a KERN ODC 251 tablet microscope camera, see table for details
- The high-resolution KERN ODC 861 20 MP fluorescence camera with cooling enables razor-sharp and low-noise images
- The intuitive KERN ODC 251 microscope tablet camera combines an Android tablet with a 10.5" HD display and a high-resolution 8 MP 4K camera
- The C-mount adapter provides the right connection between microscope and camera, enabling stable, high-resolution image transmission

Technical data

- OBN 142**
- Quintuple nosepiece
 - Infinity optical system
 - Siedentopf 30° inclined/360° rotatable
 - Diopter adjustment, both-sided
 - Eyepiece: HWF 10x/Ø 20 mm
 - Overall dimensions W×D×H 390×235×620 mm
 - Net weight approx. 11 kg
- ODC 861**
- 20 MP fluorescence camera with Peltier cooling
 - Maximum resolution 5440×3648 Pixel
 - Microscope software VIS Pro, see *internet*
- ODC 251**
- 10,5" LCD-Touchscreen
 - Screen resolution: pixels 1920×1280
 - CPU: Quad Core Cortex-A17; 1,8 GHz
 - Measuring software (pre-installed)

For detailed information on the individual components, see the relevant product description of the individual item

Model		Included camera	Camera Resolution	Camera Interface	Camera Framerate	Camera Sensor
KERN						
OBN 142C861	NEW	ODC 861	20 MP	USB 3.0	5 - 30 fps	Sony CMOS 1"
OBN 142T251	NEW	ODC 251	8 MP 4 K	Ethernet, Bluetooth, WiFi, USB 2.0, USB 3.0, HDMI	30 fps	CMOS 1/1,8"

NEW New model