


## PRODUCT BROCHURE

---

# KERN 318-06

Date: 10.03.2026

ENGLISH



Milligram weights, E2, 50 mg, Platelet, stainless steel polished

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- Milligram weights, E2, 50 mg, Platelet, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 318-06                         |
| Model-Serie           | 318                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML E2 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 404576 1240452                 |
| Predecessor Model     | 318-06-2019a                   |
| Customs tariff number | 84239090                       |
| Country of Origin     | CN                             |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 50 mg       |
| Tolerance (OIML) | 0,012 mg    |
| Case aluminium   | 317-009-600 |
| Case plastic     | 347-009-400 |
| Case wood        | 338-090-200 |

### Technical data - Construction

|                    |                          |
|--------------------|--------------------------|
| Dimensions (W×D×H) | 7,5×6×3 mm               |
| Material           | stainless steel polished |

### Technical data - Measuring system

|                              |          |
|------------------------------|----------|
| Nominal value                | 50 mg    |
| OIML class                   | E2       |
| OIML conform                 | yes      |
| Tolerance (OIML)             | 0,012 mg |
| Density (kg/m <sup>3</sup> ) | 7960     |
| Volume (cm <sup>3</sup> )    | 0,0058   |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 7.5×6×3 mm     |
| Gross weight                 | 0.05 kg        |
| Net weight                   | 0.05 kg        |
| Shipping method              | Parcel service |
| Delivery time                | 1 d            |

### Technical data - Verification

Verification possible      yes

## Optional Services

DAkkS Calibration      962-356

## FUNCTIONS

Standard



Option



## ACCESSORIES

| Model                              | Description                                    |
|------------------------------------|------------------------------------------------|
| <a href="#"><u>338-090-200</u></a> | Wooden Box for Set of Weights KERN 338-090-200 |
| <a href="#"><u>347-009-400</u></a> | Plastic Weight Box KERN 347-009-400            |
| <a href="#"><u>317-009-600</u></a> | Aluminium Protective Box KERN 317-009-600      |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242                          |
| <a href="#"><u>318-270</u></a>     | Dusting Brush KERN 318-270                     |

## SERVICES

| Model                           | Description                                          |
|---------------------------------|------------------------------------------------------|
| <a href="#"><u>962-356</u></a>  | DAkkS Calibration KERN 962-356                       |
| <a href="#"><u>962-356R</u></a> | DAkkS Recalibration KERN 962-356R                    |
| <a href="#"><u>921-356</u></a>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-356 |